



COMPUTER SPECIFICATIONS

Intel i5 or better (NO Small Form Factor (SFF) or Micro computers) * NO VIA or SIS chipsets Regular or Active USB extension cable to connect sensors to the pc.
Regular or Active USB extension cable to connect sensors to the no
Regular of Active God extension came to connect sensors to the pc.
SFF pcs are usually low powered and lack space to add additional cards if needed.
Vorkstation: Windows 10, Windows 11 – Professional versions
Vindows Server Versions: 2016 / 2019 / 2022
Check compatibility of hardware drivers with version of Windows.
Windows versions should be Microsoft Supported versions
Workstation: 8 GB or better
Server: 16 GB or better
More RAM increases performance
Workstation: 80+ GB
Server: 500+ GB for data storage
Widescreen monitor (1366 x 768 minimum resolution)
Dedicated Video Card with 2GB+ Memory (NVidia Chipset recommended)
1 Gigabit WIRED Network Card (Wireless NOT recommended)
HIGH SPEED Internet connection for support
Windows approved backup system for CADI Patients' folders. Backups should be stored off site.